Transcript: Advisory on Systematic Underestimation of Liver Stiffness

## Contents

Transcript: Advisory on Systematic Underestimation of Liver Stiffness .....	2
Transkript: Hinweis zur systematischen Unterschätzung der Lebersteifigkeit (DE) .....	4
Transcripción: Aviso sobre la subestimación sistemática de la rigidez hepática (ES) .....	6

## Transcript: Advisory on Systematic Underestimation of Liver Stiffness

### Slide 1 – Introduction

This is a short training video to summarize an advisory related to systematic underestimation of liver stiffness due to small pixel size.

### Slide 2 – Summary of Risk

When smaller-than-recommended pixel sizes—those less than 1.65 millimeters—are used with MMDI versions 3.0 or 3.1, liver stiffness values may be systematically underestimated.

This may lead to underestimation of fibrosis stage.

The error is not visually apparent.

The issue can be fully corrected by adjusting the protocol to achieve a 1.65 millimeter pixel size, or a minimum of 1.5 millimeters on Philips systems.

### Slide 3 – Description of the Problem

The multi-modal direct inversion (MMDI) algorithm is used to create standard contrasts for MRE, including stiffness, wave, and confidence.

Stiffness of a pixel is estimated using data in a small, circular kernel.

In MMDI versions 3.0 and 3.1, stiffness is calculated using a fixed 11-pixel kernel.

Because the kernel size is defined in pixels rather than millimeters, changes in image pixel size alter the physical size of the kernel.

### Slide 4 – What Does This Mean?

A 420 mm field of view with a 256x256 matrix produces 1.65 mm pixels (validated standard).

Smaller pixels (<1.65 mm) result in lower stiffness values.

Larger pixels (>1.65 mm) result in higher stiffness values.

This creates a measurement bias directly linked to pixel size.

## Slide 5 – Affected Configurations

Applies to MMDI versions 3.0 and 3.1.

Observed in Philips MR Elastography with default acquisition protocols (1.2–1.4 mm voxel size).

## Slide 6 – Clinical Impact

Bias may underestimate liver stiffness.

Patients may be classified one or two fibrosis stages lower than actual.

This can impact clinical decisions and treatment planning.

## Slide 7 – Detectability

Images appear normal and high quality.

Error is not visually detectable.

Detection requires reviewing pixel size (DICOM tag) and MMDI version.

## Slide 8 – Final Summary

Smaller pixel sizes can lead to underestimated stiffness and fibrosis stage.

Error is not visible.

Correctable by adjusting protocol to 1.65 mm pixel size.

Future MMDI 3.2 will include adaptive kernel scaling.

## Slide 9 – Closing

More information: [resoundant.com/pixel-size-advisory](https://resoundant.com/pixel-size-advisory)

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